















Locking Screw and Drill Bit

Reference Chart

	Product Image	Item Number	Thread Diameter	Screw Type	Length Range	Drill Bit	Drill Bit Number	Tap	Tap Number
T6 Hexalobe		VAR-8920L-XX	2.0 mm	Low-Profile Locking Screw	6 mm-14 mm increments of one 16 -30 mm increments of two	1.5 mm	VAR-4020D	2.0 mm	VAR-4020T
T8 Hexalobe		VAR-8924L-XX	2.4 mm	Low-Profile Locking Screw	8 mm-14 mm increments of one 16 mm-30 mm increments of two	1.8 mm	VAR-4024D	2.4 mm	VAR-4024T
T10 Hexalobe		VAR-8827L-XX	2.7 mm	Low-Profile Locking Screw	10mm-34mm increments of two	2.0 mm	VAR-8944-22	2.7 mm	VAR-4027T
		VAR-8930L-XX	3.0 mm	Low-Profile Locking Screw	8 mm-40 mm increments of two	2.3 mm	VAR-4030D	3.0 mm	VAR-4030T
T15 Hexalobe		VAR-8835L-XX	3.5 mm	Low-Profile Locking Screw	8 mm-60 mm increments of two	2.8 mm	VAR-4035D	N/A	N/A
		VAR-8840L-XX	4.0 mm	Low-Profile Locking Screw	18 mm-60 mm increments of two 65 mm	3.5 mm	VAR-4040D	N/A	N/A
T20 Hexalobe		VAR-8845L-XX	4.5 mm	Low-Profile Locking Screw	24 mm-65 mm increments of two 65mm	3.8 mm	VAR-4045D	N/A	N/A





Cortical Screw and Drill Bit

Reference Chart

	Product Image	Item Number	Thread Diameter	Screw Type	Length Range	Drill Bit	Drill Bit Number	Tap	Tap Number
T6 Hexalobe		VAR-8916-XX	1.6 mm	Cortical Screw	6 mm-14 mm increments of one 16 mm-20 mm increments of two	1.1 mm	VAR-4016D	N/A	N/A
		VAR-8920-XX	2.0 mm	Cortical Screw	6 mm-14 mm increments of one 16 mm-30 mm increments of two	1.5 mm	VAR-4020D	2.0 mm	VAR-4020T
T8 Hexalobe		VAR-8924-XX	2.4 mm	Cortical Screw	8 mm-14 mm increments of one 16 mm-30 mm increments of two	1.8 mm	VAR-4024D	2.4 mm	VAR-4024T
T10 Hexalobe		VAR-8827-XX	2.7 mm	Cortical Screw	10 mm-34 mm increments of two	2.0 mm	VAR-8944-22	2.7 mm	VAR-4027T
		VAR-8930-XX	3.0 mm	Cortical Screw	8 mm-40 mm increments of two	2.3 mm	VR-4030D	3.0 mm	VAR-4030T
T15 Hexalobe		VAR-8835-XX	3.5 mm	Cortical Screw	8 mm-60 mm increments of two 65mm	2.5 mm	VAR-8943-30	N/A	N/A
T20 Hexalobe		VAR-8845-XX	4.5 mm	Cortical Screw	24 mm-60 mm increments of two 65mm	3.0 mm	VAR-8970-30	N/A	N/A

Variable Angle Screw and Drill Bit

Reference Chart

	Product Image	Item Number	Thread Diameter	Screw Type	Length Range	Drill Bit	Drill Bit Number	Tap	Tap Number
T6 Hexalobe		VAR-8916V-XX	1.6 mm	Variable Angle Screw	6 mm-14 mm increments of one 16 mm-20 mm increments of two	1.1 mm	VAR-4016D	N/A	N/A
		VAR-8920V-XX	2.0 mm	Variable Angle Screw	6 mm-14 mm increments of one 16 mm-30 mm increments of two	1.5 mm	VAR-4020D	2.0 mm	VAR-4020T
T8 Hexalobe		VAR-8924V-XX	2.4 mm	Variable Angle Screw	8 mm-14 mm increments of one 16 mm-30 mm increments of two	1.8 mm	VAR-4024D	2.4 mm	VAR-4024T
T10 Hexalobe		VAR-8930V-XX	3.0 mm	Variable Angle Screw	8 mm-40 mm increments of two	2.3 mm	VAR-4030D	3.0 mm	VAR-4030T