



OrthoLine™ System Fracture Plate

Size Chart for Distal Femur Osteotomy Plate

	Product Image	Part Number	Length (mm)	Max Thickness Head (mm)	Max Thickness Shaft (mm)	Max Width Head (mm)	Max Width Shaft (mm)
Distal Femur Osteotomy Plate, Ti, Broad, Short, 1.6 mm		VAR-3116BDFOS-L	39.4	1.8	1.8	6.04	5.0
		VAR-3116BDFOS-R					
Distal Femur Osteotomy Plate, Ti, Broad, 1.6 mm		VAR-3116BDFO-L	43.5	2.1	1.8	6.04	5.0
		VAR-3116BDFO-R					
Distal Femur Osteotomy Plate, Ti, Short, 2.0 mm		VAR-3120DFOS-L	42.7	2.5	2.3	7.4	6.0
		VAR-3120DFOS-R					
Distal Femur Osteotomy Plate, Ti, 2.0 mm		VAR-3120DFO-L	46.9	3.1	2.3	7.4	6.0
		VAR-3120DFO-R					



OrthoLine™ System Fracture Plate

Size Chart for Distal Femur Osteotomy Plate

	Product Image	Part Number	Length (mm)	Max Thickness Head (mm)	Max Thickness Shaft (mm)	Max Width Head (mm)	Max Width Shaft (mm)
Distal Femur Osteotomy Plate, Ti, Short, 2.4 mm		VAR-3124DFOS-L	54.5	3.07	3.07	9.7	7.3
		VAR-3124DFOS-R					
Distal Femur Osteotomy Plate, Ti, 2.4 mm		VAR-3124DFO-L	57.3	3.6	3.07	9.7	7.3
		VAR-3124DFO-R					
Distal Femur Osteotomy Plate, Ti, 3.0 mm		VAR-3130DFO-L	72.7	3.5	3.6	10.5	9.3
		VAR-3130DFO-R					
Distal Femur Osteotomy Plate, SS, 3.5 mm		VAR-3035DFO-L	85.6	4.1	4.3	13.05	11.9
		VAR-3035DFO-R					