


Vet CUE Instrument Set

VAR-7020S



Pic.	Item Number	Qty.	Description
1	VAR-7004		Tamp, humerus CUE
2	VAR-7016		Tamp, ulna CUE
3	VAR-7007-L		Implant Inserter, ulna CUE, large
4	VAR-7007-M		Implant Inserter, ulna CUE, medium
5	VAR-7003-L		Implant Inserter, humerus CUE, large
6	VAR-7003-M		Implant Inserter, humerus CUE, medium
7	VAR-7006-L		Trial, ulna CUE, large
8	VAR-7006-M		Trial, ulna CUE, medium
9	VAR-7014-L		Drill, ulna CUE, large
10	VAR-7014-M		Drill, ulna CUE, medium

Pic.	Item Number	Qty.	Description
11	VAR-7015-L		Guide, ulna CUE, large
12	VAR-7015-M		Guide, ulna CUE, medium
13	VAR-7002-L		Trial, humerus CUE, large
14	VAR-7002-M		Trial, humerus CUE, medium
15	VAR-7011-L		Drill, humerus CUE, large
16	VAR-7011-M		Drill, humerus CUE, medium
17	VAR-7012-L		Drill Stop, humerus CUE, large
18	VAR-7012-M		Drill Stop, humerus CUE, medium
19	VAR-7001-L		Template, humerus CUE, large
20	VAR-7001-M		Template, humerus CUE, medium

Vet CUE Instrument Set

VAR-7020S

Date of Surgery:

Surgeon Name:

Account:

Representative Contact Information:

Set Total:

Miscellaneous Notes:

To obtain a copy of a product's directions for use (DFU), which includes sterilization instructions, please visit our website at www.arthrex.com/edfu or contact your local representative.

Page 2 of 2

www.arthrexvetsystems.com

© 2019 Arthrex, Inc. All rights reserved. DOC1-000090-en-US_B



Arthrex
Vet Systems