











Distal Femur Osteotomy Plate Size Chart









OrthoLine™ System Fracture Plate

	Product Image	Part Number	Length (mm)	Max Thickness Head (mm)	Max Thickness Shaft (mm)	Max Width Head (mm)	Max Width Shaft (mm)
Distal Femur Osteotomy Plate, Ti, Broad, Short, 1.6 mm		VAR-3116BDFOS-L	39.4	1.8	1.8	6.04	5.0
		VAR-3116BDFOS-R					
Distal Femur Osteotomy Plate, Ti, Broad, 1.6 mm		VAR-3116BDFO-L	43.5	2.1	1.8	6.04	5.0
		VAR-3116BDFO-R					
Distal Femur Osteotomy Plate, Ti, Short, 2.0 mm		VAR-3120DFOS-L	42.7	2.5	2.3	7.4	6.0
		VAR-3120DFOS-R					
Distal Femur Osteotomy Plate, Ti, 2.0 mm		VAR-3120DFO-L	46.9	3.1	2.3	7.4	6.0
		VAR-3120DFO-R					



Distal Femur Osteotomy Plate Size Chart

OrthoLine™ System Fracture Plate

	Product Image	Part Number	Length (mm)	Max Thickness Head (mm)	Max Thickness Shaft (mm)	Max Width Head (mm)	Max Width Shaft (mm)
Distal Femur Osteotomy Plate, Ti, Short, 2.4 mm		VAR-3124DFOS-L	54.5	3.07	3.07	9.7	7.3
		VAR-3124DFOS-R					
Distal Femur Osteotomy Plate, Ti, 2.4 mm		VAR-3124DFO-L	57.3	3.6	3.07	9.7	7.3
		VAR-3124DFO-R					
Distal Femur Osteotomy Plate, Ti, 3.0 mm		VAR-3130DFO-L	72.7	3.5	3.6	10.5	9.3
		VAR-3130DFO-R					
Distal Femur Osteotomy Plate, SS, 3.5 mm		VAR-3035DFO-L	85.6	4.1	4.3	13.05	11.9
		VAR-3035DFO-R					