






# OrthoLine™ Fracture Management System

## Weight Reference Chart – Distal Femoral Plate



	1 kg-5 kg	4 kg-8 kg	7 kg-12 kg	13 kg-19 kg	20 kg-32 kg
<b>Plate Size</b>	1.6B mm	1.6B mm/2.0 mm	2.4 mm	3.0 mm	3.5 mm
<b>Plate</b>					

All weight ranges are suggestions based on the ultimate load of the product and/or on our veterinary consultants' recommendations. It is up to the veterinarian to determine the proper product, application, technique, and weight range for using the selected product.