











# OrthoLine™ Fracture Management System

## Weight Reference Chart – Iliac Fracture Plate

					
	1-5 kg	4-8 kg	7-12 kg	13-19 kg	20-32 kg
Plate Size	1.6B mm	1.6B mm/2.0 mm	2.4 mm	3.0 mm	3.5 mm
Plate					

All weight ranges are suggestions based on the ultimate load of the product and/or on our veterinary consultants' recommendations. It is up to the veterinarian to determine the proper product, application, technique, and weight range for using the selected product.