



OrthoLine™ System Fracture Plate

Size Chart for Distal Femur Fracture Plate

	Product Image	Part Number	Length (mm)	Max Thickness Head (mm)	Max Thickness Shaft (mm)	Max Width Head (mm)	Max Width Shaft (mm)
Distal Femur Plate, Ti, Broad, 1.6 mm		VAR-3116BDF-L	68.3	1.9	1.8	6.1	5.0
		VAR-3116BDF-R					
Distal Femur Plate, Ti, 2.0 mm		VAR-3120DF-L	81.6	3.1	2.3	7.4	6.0
		VAR-3120DF-R					
Distal Femur Plate, Ti, 2.4 mm		VAR-3124DF-L	104.4	3.5	2.8	9.6	7.3
		VAR-3124DF-R					
Distal Femur Plate, Ti, 3.0 mm		VAR-3130DF-R	131.8	3.6	3.6	10.6	9.3
		VAR-3130DF-L					



OrthoLine™ System Fracture Plate

Size Chart for Distal Femur Fracture Plate

	Product Image	Part Number	Length (mm)	Max Thickness Head (mm)	Max Thickness Shaft (mm)	Max Width Head (mm)	Max Width Shaft (mm)
Distal Femur Plate, SS, 3.5 mm		VAR-3035DF-L	161.3	4.2	4.3	13.05	11.9
		VAR-3035DF-R					